


Product Datasheet

GRASP1 (GRIPAP1) antibody (orb1338095)

Description	Mouse monoclonal antibody to GRASP1 (GRIPAP1)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human GRIPAP1 (NP_064522.3) produced in HEK293T cell.
Target	GRIPAP1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	96.0 kDa
Uniprot ID	Q4V328
NCBI	NM_020137
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...