


Product Datasheet

ATP5F1B antibody (orb1338100)

Description	Mouse monoclonal antibody to ATP5F1B
Species/Host	Mouse
Reactivity	Canine, Human, Monkey
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human ATP5B (NP_001677) produced in HEK293T cell.
Target	ATP5F1B
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	56.6 kDa
Uniprot ID	P06576
NCBI	NM_001686
Dilution Range	WB 1:1000, IHC 1:50, IF 1:100, Flow 1:100
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
Jurkat cells...