


Product Datasheet

B Raf (BRAF) antibody (orb1338104)

Description	Mouse monoclonal antibody to B Raf (BRAF)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IP, WB
Immunogen	Full length human recombinant protein of human BRAF (NP_004324) produced in HEK293T cell.
Target	BRAF
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	84.4 kDa
Uniprot ID	P15056
NCBI	NM_004333
Dilution Range	WB 1:2000, IF 1:100, Flow 1:100
Expiration Date	12 months from date of receipt.



Flow
cytometric
HEK293T
cells
transfecte...