


Product Datasheet

Claudin 5 (CLDN5) antibody (orb1338106)

Description	Mouse monoclonal antibody to Claudin 5 (CLDN5)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF
Immunogen	Full length human recombinant protein of human CLDN5(NP_003268) produced in HEK293T cell.
Target	CLDN5
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	23 kDa
Uniprot ID	O00501
NCBI	NM_003277
Dilution Range	IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



immunofluorescent staining of COS7 cells...