


Product Datasheet

c-Jun (JUN) antibody (orb1338107)

Description	Mouse monoclonal antibody to c-Jun (JUN)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Full length human recombinant protein of human JUN (NP_002219) produced in HEK293T cell.
Target	JUN
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	35.5 kDa
Uniprot ID	P05412
NCBI	NM_002228
Dilution Range	WB 1:2000, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...