


Product Datasheet

Aconitase 2 (ACO2) antibody (orb1338108)

Description	Mouse monoclonal antibody to Aconitase 2 (ACO2)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Full length human recombinant protein of human ACO2 (NP_001089) produced in HEK293T cell.
Target	ACO2
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	85.4 kDa
Uniprot ID	Q99798
NCBI	NM_001098
Dilution Range	WB 1:2000, IF 1:100, Flow 1:100
Expiration Date	12 months from date of receipt.



immunofluorescent staining of COS7 cells...