


Product Datasheet

AKT2 antibody (orb1338118)

Description	Mouse monoclonal antibody to AKT2
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Full length human recombinant protein of human AKT2 (NP_001617) produced in HEK293T cell.
Target	AKT2
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	55.6 kDa
Uniprot ID	P31751
NCBI	NM_001626
Dilution Range	WB 1:2000, IP 2ug/500ul
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...