


Product Datasheet

DCXR antibody (orb1338119)

Description	Mouse monoclonal antibody to DCXR
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human DCXR (NP_057370) produced in HEK293T cell.
Target	DCXR
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	25.9 kDa
Uniprot ID	Q7Z4W1
NCBI	NM_016286
Dilution Range	WB 1:2000, IHC 1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...