


Product Datasheet

GBE1 antibody (orb1338120)

Description	Mouse monoclonal antibody to GBE1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, IP, WB
Immunogen	Full length human recombinant protein of human GBE1 (NP_000149) produced in HEK293T cell.
Target	GBE1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	80.3 kDa
Uniprot ID	Q04446
NCBI	NM_000158
Dilution Range	WB 1:500~2000, IHC 1:50, IF 1:50~100, FLOW 1:100, IP 2ug/500ul
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...