

Product Datasheet

HSP70-1A (HSPA1A) antibody (orb1338122)

Description Mouse monoclonal antibody to HSP70-1A (HSPA1A)

Species/Host Mouse

Reactivity Canine, Human, Monkey, Mouse, Rat

Conjugation Unconjugated

Tested Applications IHC, WB

Immunogen Full length human recombinant protein of human HSPA1A (NP_005336) produced in HEK293T cell.

Target HSPA1A

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG1

Clonality Monoclonal


MW 70.0 kDa

Uniprot ID [P0DMV8](#)

NCBI [NM_005345](#)

Dilution Range WB 1:2000, IHC 1:50

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...