


Product Datasheet

PASK antibody (orb1338123)

Description	Mouse monoclonal antibody to PASK
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human PASK (NP_055963) produced in HEK293T cell.
Target	PASK
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	142.7 kDa
Uniprot ID	Q96RG2
NCBI	NM_015148
Dilution Range	WB 1:2000, IHC 1:50, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...