

## Product Datasheet

### Aconitase 2 (ACO2) antibody (orb1338124)

**Description** Mouse monoclonal antibody to Aconitase 2 (ACO2)

**Species/Host** Mouse

**Reactivity** Human, Mouse, Rat

**Conjugation** Unconjugated

**Tested Applications** IF, IHC, WB

**Immunogen** Full length human recombinant protein of human ACO2 (NP\_001089) produced in HEK293T cell.

**Target** ACO2

**Note** For research use only

**Application notes** For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

**Isotype** IgG1

**Clonality** Monoclonal


**MW** 85.4 kDa

**Uniprot ID** [Q99798](#)

**NCBI** [NM\\_001098](#)

**Dilution Range** WB 1:2000, IHC 1:50, IF 1:100

**Expiration Date** 12 months from date of receipt.



Western blot analysis of HEK293T cells w...