


Product Datasheet

Prolactin (PRL) antibody (orb1338131)

Description	Mouse monoclonal antibody to Prolactin (PRL)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Full length human recombinant protein of human PRL (NP_000939) produced in HEK293T cell.
Target	PRL
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	25.9 kDa
Uniprot ID	P01236
NCBI	NM_000948
Dilution Range	WB 1:2000, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...