


Product Datasheet

L1CAM antibody (orb1338148)

Description	Mouse monoclonal antibody to L1CAM
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, IP, WB
Immunogen	Full length human recombinant protein of human L1CAM (NP_000416) produced in HEK293T cell.
Target	L1CAM
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	140.0 kDa
Uniprot ID	P32004
NCBI	NM_000425
Dilution Range	WB 1:1000~2000, IHC 1:50, IF 1:100, FLOW 1:100, IP 2-4ug/mg
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...