


Product Datasheet

TTLL12 antibody (orb1338158)

Description	Mouse monoclonal antibody to TTLL12
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human TTLL12 (NP_055955) produced in HEK293T cell.
Target	TTLL12
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	74.4 kDa
Uniprot ID	Q14166
NCBI	NM_015140
Dilution Range	WB 1:2000, IHC 1:50, IF 1:100, Flow 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...