


Product Datasheet

HDAC10 antibody (orb1338160)

Description	Mouse monoclonal antibody to HDAC10
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IP, WB
Immunogen	Full length human recombinant protein of human HDAC10 (NP_114408) produced in HEK293T cell.
Target	HDAC10
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	71.4 kDa
Uniprot ID	Q969S8
NCBI	NM_032019
Dilution Range	WB 1:2000, Flow 1:100, IP 2-4ug/mg
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...