


Product Datasheet

RPA34 (RPA2) antibody (orb1338161)

Description	Mouse monoclonal antibody to RPA34 (RPA2)
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, IP, WB
Immunogen	Full length human recombinant protein of human RPA2 (NP_002937) produced in HEK293T cell.
Target	RPA2
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	29.1 kDa
Uniprot ID	P15927
NCBI	NM_002946
Dilution Range	WB 1:1000~2000, IHC 1:50, IF 1:100, FLOW 1:100, IP 2ug/500ul
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...