


Product Datasheet

IBRDC2 (RNF144B) antibody (orb1338231)

Description	Mouse monoclonal antibody to IBRDC2 (RNF144B)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-256 of human RNF144B (NP_877434) produced in E. coli.
Target	RNF144B
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	33.5 kDa
Uniprot ID	Q7Z419
NCBI	NM_182757
Dilution Range	IF 1:50~100, FLOW 1:100
Expiration Date	12 months from date of receipt.



immunofluorescent staining of COS7 cells...