


Product Datasheet

Sorbitol Dehydrogenase (SORD) antibody (orb1338239)

Description	Mouse monoclonal antibody to Sorbitol Dehydrogenase (SORD)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Full length human recombinant protein of human SORD (NP_003095) produced in HEK293T cell.
Target	SORD
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	38.3 kDa
Uniprot ID	Q00796
NCBI	NM_003104
Dilution Range	WB 1:1000, IF 1:100, Flow 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...