

Product Datasheet

Protein Kinase A regulatory subunit I alpha (PRKAR1A) antibody (orb1338292)

Description Mouse monoclonal antibody to Protein Kinase A regulatory subunit I alpha (PRKAR1A)

Species/Host Mouse

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications IF, IHC, WB

Immunogen Full-length protein expressed in 293T cell transfected with human PRKAR1A expression vector

Target PRKAR1A

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG1

Clonality Monoclonal


MW 43.0 kDa

Uniprot ID [P10644](#)

NCBI [NM_212472](#)

Dilution Range WB 1:2000, IHC 1:50, IF 1:100

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...