


Product Datasheet

PPP1R7 antibody (orb1338317)

Description	Mouse monoclonal antibody to PPP1R7
Species/Host	Mouse
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human PPP1R7 (NP_002703) produced in HEK293T cell.
Target	PPP1R7
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	41.4 kDa
Uniprot ID	Q15435
NCBI	NM_002712
Dilution Range	WB 1:2000, IHC 1:50, IF 1:50~100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...