

Product Datasheet

Cytochrome P450 Reductase (POR) antibody (orb1338329)

Description Mouse monoclonal antibody to Cytochrome P450 Reductase (POR)

Species/Host Mouse

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications FC, IF, IHC, WB

Immunogen Full length human recombinant protein of human POR (NP_000932) produced in HEK293T cell.

Target POR

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG1

Clonality Monoclonal


MW 76.9 kDa

Uniprot ID [P16435](#)

NCBI [NM_000941](#)

Dilution Range WB 1:2000, IHC 1:50, IF 1:50~100, FLOW 1:100

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...