


Product Datasheet

PP5 (PPP5C) antibody (orb1338345)

Description	Mouse monoclonal antibody to PP5 (PPP5C)
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human PPP5C (NP_006238) produced in HEK293T cell.
Target	PPP5C
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	56.9 kDa
Uniprot ID	P53041
NCBI	NM_006247
Dilution Range	WB 1:1000~2000, IHC 1:50, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...