


Product Datasheet

DHFR antibody (orb1338387)

Description	Mouse monoclonal antibody to DHFR
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full-length protein expressed in 293T cell transfected with human DHFR expression vector
Target	DHFR
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	21.5 kDa
Uniprot ID	P00374
NCBI	NM_000791
Dilution Range	WB 1:2000 IHC 1:50 IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...