


Product Datasheet

SCY1 like 3 (SCYL3) antibody (orb1338440)

Description	Mouse monoclonal antibody to SCY1 like 3 (SCYL3)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full-length protein expressed in 293T cell transfected with human SCYL3 expression vector
Target	SCYL3
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	82.8 kDa
Uniprot ID	Q8IZE3
NCBI	NM_020423
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...