


Product Datasheet

CD4 antibody (orb1338452)

Description	Mouse monoclonal antibody to CD4
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human CD4 (NP_000607) produced in HEK293T cell.
Target	CD4
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	51.1 kDa
Uniprot ID	P01730
NCBI	NM_000616
Dilution Range	WB 1:1000~2000, IHC 1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...