


Product Datasheet

PARP16 antibody (orb1338473)

Description	Mouse monoclonal antibody to PARP16
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Full-length protein expressed in 293T cell transfected with human PARP16 expression vector .
Target	PARP16
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	36.4 kDa
Uniprot ID	Q8N5Y8
NCBI	NM_017851
Dilution Range	FLOW 1:100
Expiration Date	12 months from date of receipt.



Flow
cytometric
Flow
cytometric
HEK293T ...