

Product Datasheet

Melanoma gp100 (PMEL) antibody (orb1338520)

Description Mouse monoclonal antibody to Melanoma gp100 (PMEL)

Species/Host Mouse

Reactivity Canine, Human, Monkey, Rat

Conjugation Unconjugated

Tested Applications IF, IHC, WB

Immunogen Full-length protein expressed in 293T cell transfected with human SILV expression vector

Target PMEL

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG2b

Clonality Monoclonal


MW 70.3 kDa

Uniprot ID [P40967](#)

NCBI [NM_006928](#)

Dilution Range WB 1:1000, IHC: 1:150

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...