


Product Datasheet

PLK1 antibody (orb1338527)

Description	Mouse monoclonal antibody to PLK1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC, IP, WB
Immunogen	Full length human recombinant protein of human PLK1 (NP_005021) produced in HEK293T cell.
Target	PLK1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	68.1 kDa
Uniprot ID	P53350
NCBI	NM_005030
Dilution Range	WB: 1:200 - 1:1000, IHC 1:150, IP 2ug/500ul
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...