


Product Datasheet

TrkB (NTRK2) antibody (orb1338535)

Description	Mouse monoclonal antibody to TrkB (NTRK2)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full-length protein expressed in 293T cell transfected with human NTRK2 expression vector
Target	NTRK2
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	92 kDa
Uniprot ID	Q16620
NCBI	NM_001007097
Dilution Range	WB 1:1000, IHC 1:150, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...