


Product Datasheet

Chk2 (CHEK2) antibody (orb1338537)

Description	Mouse monoclonal antibody to Chk2 (CHEK2)
Species/Host	Mouse
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC, IP, WB
Immunogen	Full length human recombinant protein of human CHEK2 (NP_009125) produced in HEK293T cell.
Target	CHEK2
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	60.7 kDa
Uniprot ID	O96017
NCBI	NM_007194
Dilution Range	WB 1:1000~2000, IHC 1:150, FLOW 1:100, IP 2ug/500ul
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...