


Product Datasheet

SIGLEC9 antibody (orb1338539)

Description	Mouse monoclonal antibody to SIGLEC9
Species/Host	Mouse
Reactivity	Canine, Human, Monkey
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human SIGLEC9 (NP_055256) produced in HEK293T cell.
Target	SIGLEC9
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	50.1 kDa
Uniprot ID	Q9Y336
NCBI	NM_014441
Dilution Range	WB 1:500~2000, IHC 1:50, IF 1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...