


Product Datasheet

TRPM4 antibody (orb1338540)

Description	Mouse monoclonal antibody to TRPM4
Species/Host	Mouse
Reactivity	Canine, Human, Monkey
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human TRPM4 (NP_060106) produced in HEK293T cell.
Target	TRPM4
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	134.3 kDa
Uniprot ID	Q8TD43
NCBI	NM_017636
Dilution Range	WB 1:5000~10000, IHC 1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...