


Product Datasheet

TRPM4 antibody (orb1338544)

Description	Mouse monoclonal antibody to TRPM4
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full-length protein expressed in 293T cell transfected with human TRPM4 expression vector
Target	TRPM4
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	134.3 kDa
Uniprot ID	Q8TD43
NCBI	NM_017636
Dilution Range	WB 1:1000, IHC 1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...