


Product Datasheet

alpha 1 Antitrypsin (SERPINA1) antibody (orb1338546)

Description	Mouse monoclonal antibody to alpha 1 Antitrypsin (SERPINA1)
Species/Host	Mouse
Reactivity	Canine, Human, Monkey
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	Full length human recombinant protein of human SERPINA1 (NP_000286) produced in HEK293T cell.
Target	SERPINA1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	44.3 kDa
Uniprot ID	P01009
NCBI	NM_000295
Dilution Range	WB 1:500~1000, IHC 1:150, IP 2ug/500ul
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...