

Product Datasheet

alpha 1 Antitrypsin (SERPINA1) antibody (orb1338547)

Description Mouse monoclonal antibody to alpha 1 Antitrypsin (SERPINA1)

Species/Host Mouse

Reactivity Canine, Human, Monkey

Conjugation Unconjugated

Tested Applications FC, IHC, IP, WB

Immunogen Full-length protein expressed in 293T cell transfected with human SERPINA1 expression vector

Target SERPINA1

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG1

Clonality Monoclonal


MW 46.7 kDa

Uniprot ID [P01009](#)

NCBI [NM_000295](#)

Dilution Range WB 1:500~1000, IHC 1:150, FLOW 1:100, IP 2ug/500ul

Expiration Date 12 months from date of receipt.



Flow
cytometric
Flow
cytometric
HEK293T ...