


Product Datasheet

Prominin 2 (PROM2) antibody (orb1338564)

Description	Mouse monoclonal antibody to Prominin 2 (PROM2)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Full-length protein expressed in 293T cell transfected with human PROM2 expression vector
Target	PROM2
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	91.883 kDa
Uniprot ID	Q8N271
NCBI	NM_144707
Dilution Range	WB 1:2000, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...