


Product Datasheet

Neurogenin 1 (NEUROG1) antibody (orb1338565)

Description	Mouse monoclonal antibody to Neurogenin 1 (NEUROG1)
Species/Host	Mouse
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IP, WB
Immunogen	Full length human recombinant protein of human NEUROG1 (NP_006152) produced in HEK293T cell.
Target	NEUROG1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	25.7 kDa
Uniprot ID	Q92886
NCBI	NM_006161
Dilution Range	WB 1:500~1000, IF 1:50, IP 2-4ug/mg
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...