


Product Datasheet

Neurogenin3 (NEUROG3) antibody (orb1338566)

Description	Mouse monoclonal antibody to Neurogenin3 (NEUROG3)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, IP, WB
Immunogen	Full-length protein expressed in 293T cell transfected with human NEUROG3 expression vector
Target	NEUROG3
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	23.077 kDa
Uniprot ID	Q9Y4Z2
NCBI	NM_020999
Dilution Range	WB 1:2000, IHC 1:50 - 1:150, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...