


Product Datasheet

GFAP antibody (orb1338571)

Description	Mouse monoclonal antibody to GFAP
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full-length protein expressed in 293T cell transfected with human GFAP expression vector
Target	GFAP
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	49.9 kDa
Uniprot ID	P14136
NCBI	NM_002055
Dilution Range	WB: 1:200 - 1:1000, IHC 1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...