


Product Datasheet

GFAP antibody (orb1338574)

Description	Mouse monoclonal antibody to GFAP
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Full length human recombinant protein of human GFAP (NP_002046) produced in HEK293T cell.
Target	GFAP
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	49.9 kDa
Uniprot ID	P14136
NCBI	NM_002055
Dilution Range	WB 1:500~1000, IHC 1:50, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...