


Product Datasheet

Ki67 (MKI67) antibody (orb1338603)

Description	Mouse monoclonal antibody to Ki67 (MKI67)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant protein expressed in E.coli corresponding to amino acids 1-250 of human MKI67.
Target	MKI67
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	28.9 kDa
Uniprot ID	P46013
NCBI	NM_002417
Dilution Range	WB 1:2000 IHC 1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...