


Product Datasheet

CDX2 antibody (orb1338608)

Description	Mouse monoclonal antibody to CDX2
Species/Host	Mouse
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant protein expressed in E.coli corresponding to amino acids 1-313 of human CDX2
Target	CDX2
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	33.5 kDa
Uniprot ID	Q99626
NCBI	NM_001265
Dilution Range	WB 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...