


Product Datasheet

Oncostatin M (OSM) antibody (orb1338623)

Description	Mouse monoclonal antibody to Oncostatin M (OSM)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant protein expressed in E.coli corresponding to amino acids 25-252 of human OSM
Target	OSM
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	25.7 kDa
Uniprot ID	P13725
NCBI	NM_020530
Dilution Range	WB 1:1000, IHC 1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...