


Product Datasheet

Carrier-free (BSA/glycerol-free) DDK (FLAG) antibody (orb1338638)

Description	Mouse monoclonal antibody to Carrier-free (BSA/glycerol-free) DDK (FLAG)
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, IP, WB
Immunogen	Anti-DDK monoclonal antibody is produced by immunizing mice with a synthetic peptide coupled to KLH
Target	DDK (FLAG)
Preservatives	Lyophilized powder (original buffer 1X PBS, pH 7.3, 8% trehalose)
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
Dilution Range	WB: 1:2000, IP: 1:200, IF 1:1000, FC: 1:1000, IHC: 1:1000-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of Immunoblot anal...