


Product Datasheet

Apolipoprotein M (APOM) antibody (orb1338651)

Description	Mouse monoclonal antibody to Apolipoprotein M (APOM)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant protein expressed in E.coli corresponding to amino acids 23-188 of human APOM
Target	APOM
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	21.1 kDa
Uniprot ID	O95445
NCBI	NM_019101
Dilution Range	WB: 1:1000, IHC 1;50
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...