


Product Datasheet

ICAD (DFFA) antibody (orb1338662)

Description	Mouse monoclonal antibody to ICAD (DFFA)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant protein expressed in E.coli corresponding to amino acids 1-331 of human DFF45
Target	DFFA
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	36.3 kDa
Uniprot ID	O00273
NCBI	NM_004401
Dilution Range	WB: 1:500, IF (1:100)
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...