


Product Datasheet

FOXA1 antibody (orb1338668)

Description	Mouse monoclonal antibody to FOXA1
Species/Host	Mouse
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human FOXA1 (NP_004487) produced in E.coli.
Target	FOXA1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	49.0 kDa
Uniprot ID	P55317
NCBI	NM_004496
Dilution Range	WB 1:500~1000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...