


Product Datasheet

Neurotrophin 4 (NTF4) antibody (orb1338685)

Description	Mouse monoclonal antibody to Neurotrophin 4 (NTF4)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 81-210 human NT4 (NP_006170) produced in E.coli.
Target	NTF4
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	22.2 kDa
Uniprot ID	P34130
NCBI	NM_006179
Dilution Range	WB 1:20000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...