


Product Datasheet

Cytokeratin 19 (KRT19) antibody (orb1338698)

Description	Mouse monoclonal antibody to Cytokeratin 19 (KRT19)
Species/Host	Mouse
Reactivity	Canine, Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 240-390 of human CK19 (NP_002267) produced in E.coli.
Target	KRT19
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	43.9 kDa
Uniprot ID	P08727
NCBI	NM_002276
Dilution Range	WB 1:1000~2500, IHC 1:50, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...